iMicro Q
User Guide

Components:

iMicro Q

Shim pad

Reflecting pad

μ Ruler

Storage case/stand

Introduction

iMicro Q works with your smartphone to observe tiny objects and record digital micrographs as a portable digital microscope.

Install the APP

Scan the QR code to download and install the APP i-SEEING. For Android version, you need to set the permission by going to Setting->permission setting->App permission setting->i-Seeing to set camera to be accessible.

Shanghai QingYing E&T LLC

iOS App at Apple App Store


Android APP at Google Play

Shanghai QingYing E&T LLC

Apps are also accessible from website (www.mobile-scope.com).

Install iMicro Q

Pull it from the storage case, align the lens with phone camera and the light guide with flash, attach it to the phone through its nano-suction pad.

For phones with extruded camera, shim pad can compensate the height difference between camera and phone surface.

Make observation


Turn on the flash light, put iMicro Q against the objects to observe, the phone will make autofocus and clear image will appear on screen.

The storage case can be used as a stand to assist observation.

When observe transparent samples, such as μ Ruler or sample slide, the reflecting pad can be


Shanghai QingYing E&T LLC

used to reflect light to illuminate the sample.

The μ Ruler

5 types of ruler can help you in the calibration of the scale bar in APP, to enhance the capability with size measurement.

- A. 0.1mm line ruler
- B. 0.05mm line ruler
- C. 0.1mm cross ruler
- D. 0.5mm grid ruler
- E. 0.2mm grid ruler

Calibrate the scale bar

If you want to measure the size of micro-objects, you need to calibrate the scale bar.

Make observation of the μ Ruler, enter the calibration mode* by pushing the gear button

*The calibration code is on the label on the storage case, under the bar code.

and tune the magnification number to match the scale bar with the real length of the μ Ruler.

Maintenance

- ◆ Clean the surface with soft cloth, no organic solvent should be used.
- ◆ Keep the lens free of dust, clean the lens with soft lens tissue.
- ◆ The nano-suction pad on the backside of the device can be cleaned with PSA tape if it is contaminated by dust.

Specification

- ◆ Magnification: ~100-800X
- ◆ Max view field: ~1mm
- ◆ Size: ~26mmX12mmX3.5mm
- ◆ Weight: ~1.5g

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©2018 Shanghai QingYing E&T LLC. All rights reserved.