iMicro Q2

iMicro C

User Guide

Components:

iMicro Q2

μ Ruler M (optional)

iMicro C

Focus ring

Reflecting pad

Storage case/stand

Introduction

iMicro Q2 (iMicro C) works with your smartphone to observe tiny objects and record digital micrographs as a portable digital microscope.

Install the APP

Scan the QR code to download and install the APP i-SEEING. For the Android version, you need to set the permission by going to Setting->permission setting->App permission setting->i-Seeing to set the camera to be accessible.

Shanghai QingYing E&T LLC

iOS App at Apple App Store

Android APP at Google Play

Shanghai QingYing E&T LLC

Apps are also accessible from the website (www.mobile-scope.com).

Install iMicro

Pull it from the storage case, align the lens with the phone's primary camera, attach it to the phone through its nano-suction pad.

Make observation


Turn on the flashlight and put iMicro Q2 against the objects to observe; the phone will make autofocus and clear image appear on the screen. The storage case can be used to reflect the flashlight.

For iMicro C, the focus ring could be attached to the ring and make a stable distance between the sample and the lens, thus for a stable and easy observation.

When observing transparent samples, such as μ Ruler or sample


Shanghai QingYing E&T LLC

slide, the reflecting pad can reflect light to illuminate the sample. The storage case can be used as a stand to assist observation.

The μ Ruler M


Five types of rulers can help you calibrate the scale bar in APP, to enhance the capability with size measurement.

- A. 0.1mm line ruler
- B. 0.05mm line ruler
- C. 0.1mm cross ruler
- D. 0.5mm grid ruler
- E. 0.2mm grid ruler

Calibrate the scale bar

If you want to measure the size of micro-objects, you need to calibrate the scale bar.

Observe the μ Ruler, enter the calibration mode* by pushing the gear button

*The calibration code is on the label on

the storage case, under the bar code.
and tune the magnification number to match the scale bar with the real length on the μ Ruler M.

Maintenance

- ◆ Clean the surface with soft cloth, no organic solvent should be used.
- ◆ Keep the lens free of dust, clean the lens with soft lens tissue.
- ◆ Please keep the iMicro and phone camera clear of dust and oil to ensure stable

attachment of the lens to phone camera.

- ♦ The nano-suction pad on the backside of the device can be cleaned with stationary tape if it is contaminated by dust.

Specification

iMicro Q2

- ◆ Magnification: ~100-800X
- ◆ Max view field: ~1.8mm
- ◆ Working distance: 0.2mm
- ◆ Size: ~12.7mmX3.2mm
- ◆ Weight: ~0.37g

iMicro C

- ◆ Magnification: ~25-200X
- ◆ Max view field: ~6mm
- ◆ Working distance: ~4mm
- ◆ Size: ~12.7mmX3.2mm
- ◆ Weight: ~0.34g

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©2021 Shanghai QingYing E&T LLC. All rights reserved.